





# iAED-M1

Automated External Defibrillator (AED)



# Product Appearance

---



# Product Advantages


---

## Small, Lightweight and Portable

- Size: 17 cm \* 16 cm (6.7 in \* 6.3 in), or about the size of a table tennis racket
- Thickness: only 4.7 cm (less than 2 in)
- Weight: only 1.1 kg, or about the weight of two bottles of water

## Easy-to-Use and Quick Rescue

- Audio instructions throughout the operation, Simple three-step operation
- Simultaneously initiate heart rhythm analysis and pre-charge, Quick rescue
- Child/Adult mode switch with a single button



## Accurate Heart Rhythm Identification, Safe and Efficient

- Low-energy defibrillation technology for less myocardial damage

Jousing adopts the low-energy biphasic truncated exponential (BTE) waveform defibrillation technology and enables the device to automatically adjust the shock parameters based on the patient's transthoracic impedance, leading to good defibrillation results with less myocardial damage

- Detection and removal of pacemaker artifacts

Able to identify a sudden cardiac arrest (SCA) in patients with a pacemaker

## Excellent Performance, Long Endurance

- $\geq 5$  years of standby time
- 10 years of main unit service life
- Resistant to water and dust, anti-fall and endurable under pressure



## Application Scenarios



Vehicles



Home



Sports



Journey



## Product Features

|                        |                           |
|------------------------|---------------------------|
| Weight                 | 1.1±0.2kg                 |
| Size                   | 170x160x47(mm)            |
| Ingress Protection     | IP55                      |
| Anti-drop              | ≥1.2m                     |
| Main Unit Life         | 10 years                  |
| Battery Life           | ≥5 years                  |
| Pads Life              | 2.5 years                 |
| Storage Temperature    | 0℃~50℃                    |
| Altitude               | -390m~5920m               |
| Wireless Communication | Bluetooth, Wi-Fi or 4G/5G |


## Accessories



Main Unit



Battery




Pads



Carrying Case

# Cases of Emergency Treatment

Till now, Jousing AEDs have helped save more than 100 SCA patients



\* Some cases are not included for privacy reasons





# A Handy Device to Guard Family

Jousing Medical Co., Ltd.

301&401, Building 21, 200 Xingpu Road, Suzhou, China

400-820-9952

[www.jousing.com](http://www.jousing.com)



AED First-aid  
Map Applet



Follow Jousing  
on WeChat